

GOLEM

Tables Couches Armchairs Chairs



EUROPEAN UNION
European Regional Development Fund
Operational Programme Enterprise and Innovations
for Competitiveness

CE
ISO 9001
Since 1991

RQL

GOLEM



Armchair is the optimal solution for long-lasting stay of patient during transfusion or dialysis. Comfortable armchair for all patient positions. Travelling version with 75 mm braked wheels. Modern design with soft upholstery and white washable plastic covers. The three-segment work area with overall dimensions 640 × 2 200 mm is covered with artificial leather in color of your choice. Working area is divided into backrest, seat and leg segment. Patient comfort will also provide a headrest whose position can be adjust for patient needed. Artificial leather is washable and disinfectable with a certificate for use in healthcare.

Standard equipment is both side use infusion stand, height adjustable foot support, holder of paper rolls, a basket for personal things of a patient and a handle for a more comfortable transport.

The armchair has removable washable infusion armrests that can be adjusted vertically in all directions and inclinations. Special mechanism ensures that the selected position of the support will be retained even when changed inclination of the backrest.

The installation dimensions of the armchair are approx. 990 × 1 500 mm in the sitting position, 990 × 2 200 mm in the lying position (including armrests). Clearance profile without armrests is 750 mm.

Load weight capacity 150 kg. Weight of chair 80 - 110 kg.

If the chair is equipped with electric drives, all positions can be simply adjust even with fully load by patient or staff. If necessary you can also adjust position of Trendelenburg. Silent electric drives work with a safe voltage of 24V. The electronics ensure their reliability and smooth running.

The power of the chair is from standard socket 230V. For extra charge it can be equipped by built-in fully automatic charger. In case of power failure, the chair will be automatically switched to built-in rechargeable batteries.

For an extra charge it can be equipped plastic washable working table that can be used from both sides or flip away of the chair for facilitating the entry and exit of the patient. It can be also equipped by LED light on flexible neck for more comfort of patient (for extra charge).



Transfusion and dialysis chair



GOLEM DIA P

Code Nr. G 09 01

Fixed height 560 mm, back inclination adjustable by gas spring -13° to $+75^{\circ}$, foot segment adjustable by gas spring $+5^{\circ}$ to -40° , fixed slope of seat $+5^{\circ}$.

GOLEM DIA

Code Nr. G 09 02

Fixed height 560 mm, back inclination adjustable by electric drive -13° to $+75^{\circ}$, foot segment adjustable by electric drive $+20^{\circ}$ to -40° , seat inclination (Trendelenburg) adjustable by electric drive 0° to 20° , drive control by hand pilot.

GOLEM DIA E

Code Nr. G 09 03

Height adjustable by electric drive 560 - 760 mm, adjustable back inclination electric drive -13° to $+75^{\circ}$, electrically adjustable foot segment $+20^{\circ}$ to -40° , Trendelenburg position adjustable by electric drive 0° to 20° , drive control by hand pilot.



Accessories:

Working table

Code Nr. D 09 02

Built-in battery with a fully automatic charger

Code Nr. D 06 40

LED lamp on flexible neck

Code Nr. D 06 05

