

ACCESSORIES
MOBILES SUPPORTS
HOSPITAL SCREENS





Versatile, adapted to different hospital units

Adaptable, many options available

Two models Wall-mounted and two mobile bases

Easy to use, thanks to its mobility



HOSPITAL SCREENS

ACCESSORIES

Rigid panel screens.

- M1 white PVC panels
- From 5 to 13 panels
- Wall mounting or two mobile bases
- Finishing profiles available in different colours: Blue, Grey, Orange and White.
- Strap for closing the screen


SCREENS DIMENSIONS

		NUMBER OF PANELS								
		5	6	7	8	9	10	11	12	13
HEIGHT	1,65 M									
LENGTH (meter)		1,365	1,575	1,785	1,995	2,205	2,415	2,625	2,835	3,045


FINISHING PROFILES WHITE - BLUE - GREY - ORANGE



REFERENCE	DESCRIPTION
TENFC	Central lock
TENENT1	Brace 10 cm
TENENT2	Brace 20 cm
TENENT3	Brace 30 cm
TENPFM	Wall fixing



Central lock
TENFC



Brace
TENENT1 2 3



Wall fixing
TENPFM

